

Allen® PAL™ Paediatric Stirrups

Allen® Junior PAL™ Stirrups, and Allen® Little PAL™ Stirrups are specifically designed for paediatric patients of different ages.

The Junior PAL and Little PAL stirrups enable surgical personnel to lift and position a paediatric patient's legs and then securely hold the legs in position while surgical procedures are performed. The Allen PAL paediatric stirrup family is unique because it contains Little PAL stirrups designed for toddlers, and Junior PAL stirrups designed for preadolescent children. Other paediatric stirrups currently on the market are only available in a single size. In addition, the stirrups are designed to reduce pressure on the patient's calf, while protecting the ankle and heel.



Ideal For Use With...

All paediatric urological and laproscopic procedures including:

- Rectal Biopsy
- Cystoscopy
- Urethral Dilation
- Endoscopic Deflux Injection
- Pull Through Surgery

About The System...

The Allen Paediatric PAL Stirrups provide secure leg positioning for paediatric patients of all ages. Simultaneous abduction and lithotomy adjustments are easily made via single handle, while maintaining the sterile field. The boot is adjusted with a single clamp, making positioning simple. They attach to the rail of any operating table.

Allen Junior PAL stirrups are appropriate for most paediatric patients ages 7-11 years old. They support patients that weigh up to 113 kg (250 lbs.). Allen Little PAL stirrups are appropriate for most paediatric patients ages 3-6 years old. They support patients that weigh up to 73 kg (160 lbs.).

Key Features & Benefits

- Boot height and width designed to accommodate paediatric patients of various sizes
- Floating boot design minimizes pressure on the calf when the stirrup is in motion
- Boot's deep heel design ensures proper weight distribution by cradling the paediatric calf, ankle, and heel in a natural position.
- Boot padding extends beyond edge of boot (at calf) to minimize pressure
- Reusable Clamshell Pad – Completely encapsulates the foot, ankle, and calf
- Visual Indicators – Help ensure precise positioning
- Quick-release Handle: Enables easy intraoperative adjustments, and secures patient into desired position.

Specifications

Product Specifications:

- Attaches to any operating theatre table
- Lithotomy Range – Set the stirrup in any position between +90° to –35° lithotomy.
- Abduction Range – Set the stirrup in any position between +25° to –9° abduction.

Allen Junior PAL Stirrups

- 113 kg patient (250 lb.) maximum patient weight

Allen Little PAL Stirrups

- 227 kg patient (160 lb.) maximum patient weight

A-10055 Allen Paediatric Stirrup Package (Includes 1 Pair Junior PAL + 1 Pair Little PAL)

A-10056 Allen Junior PAL Stirrups (1 Pair)

A-10057 Allen Little PAL Stirrups (1 Pair)

(Specify) Allen Easy Lock™ Blade Clamp* (Select Rail Type)
A-40041 (EU), A-40042 (UK), A-40043 (JPN), A-40040 (US)

(Specify) Allen Rail Clamp* (Select Rail Type)
O-RCXUK (EU/UK), O-RCXD (Denyer), O-RC29X32MM (JPN), O-RCS (Swiss), O-RC2 (US)

*Allen Paediatric stirrups can be used with the Easy Lock Blade Clamp or Allen Rail Clamp. Each stirrup requires 1 clamp.



Allen Junior PAL Stirrups



Allen Little PAL Stirrups



Allen Stirrup Cart



Allen Clamshell Pad



Allen Easy Lock Clamp



Allen Rail Clamp